




**Standard Sets and  
Sets in Graphic Cases for Trauma Surgery**



## Standard Instrument Set for Dynamic and Condylar Screws and Plates in Graphic Case

**42-T-0233**

Standard Instrument Set for Dynamic Hip and Condylar Screws and Plates



**42-T-0234**

Tray (empty) for Instrument Set, Dynamic Hip and Condylar Screws and Plates (Steel)

**Listing:**

- 42-T-0233** Standard Instrument Set for Dynamic Hip and Condylar Screws and Plates
- 42-T-0234** Tray for Instrument Set, Dynamic Hip and Condylar Screws and Plates (Steel)

<u>Cat. No.:</u>		<u>Pcs.:</u>
<b>42-02-0300</b>	Calibrated Threaded Guide Pin Dia. 2.5mm x 230 mm Length, Threaded Length 10 mm	10
<b>42-02-0302</b>	Drill Guide adjustable to 125° / 130° / 135° / 140° / 145° / 150°	1
<b>42-02-0304</b>	Measuring Gauge	1
<b>42-02-0306</b>	Triple Reamer complete for Dynamic Hip Plate System	1
<b>42-02-0308</b>	Tap	1
<b>42-02-0310</b>	Guide Sleeve, short, for tapping	1
<b>42-02-0312</b>	Guide Sleeve, long, for screwing in of Lag Screws	1
<b>42-02-0314</b>	Guide Shaft for Insertion of Coupling Screws	1
<b>42-02-0316</b>	Coupling Screw for Femoral Dynamic Hip Screws	1
<b>42-02-0317</b>	Guide Shaft	1
<b>42-02-0318</b>	Wrench	1
<b>42-02-0320</b>	Impactor	1
<b>42-02-0322</b>	Handle	1
<b>42-02-0486</b>	Drill Guide 95°	1
<b>42-02-0488</b>	Triple Reamer complete for Dynamic Condylar System	1





**70-7110-10**  
Sterilization Container (recommended)

## Dynamic Hip Screw Standard Implant Set in Graphic Case

### 42-T-2400

Dynamic Hip Screw Standard Implant Set (Steel Plates in Mat. DIN ISO 5832-1)



### 42-T-2410

Tray (empty) for Dynamic Hip Screw Standard Implant Set (Steel)

#### Listing:

**42-T-2400** Dynamic Hip Screw Standard Implant Set (Steel Plates in Mat. DIN ISO 5832-1)

**42-T-2410** Tray for Dynamic Hip Screw Standard Implant Set (Steel)

Cat. No.:	Dynamic Lag Screw, Thread Length 22 mm	Pcs.:
42-02-0406	Length 65 mm, Steel	1
42-02-0408	Length 70 mm, Steel	1
42-02-0410	Length 75 mm, Steel	1
42-02-0412	Length 80 mm, Steel	1
42-02-0414	Length 85 mm, Steel	1
42-02-0416	Length 90 mm, Steel	1
42-02-0418	Length 95 mm, Steel	1
42-02-0420	Length 100 mm, Steel	1
42-02-0422	Length 105 mm, Steel	1
42-02-0424	Length 110 mm, Steel	1
42-02-0426	Length 115 mm, Steel	1

42-02-0440	Compression Screw	3
------------	-------------------	---

Cat. No.:	Dynamic Hip Screw Plates 135°	Pcs.:
42-02-0352	4 Holes, Length 78 mm, Steel	3
42-02-0354	5 Holes, Length 94 mm, Steel	2
42-02-0356	6 Holes, Length 110 mm, Steel	2

Cat. No.:	Dynamic Hip Screw Plates 150°	Pcs.:
42-02-0382	4 Holes, Length 78 mm, Steel	1
42-02-0384	5 Holes, Length 94 mm, Steel	1
42-02-0386	6 Holes, Length 110 mm, Steel	1



**70-7110-10**  
Sterilization Container (recommended)



## Dynamic Condylar Standard Implant Set in Graphic Case

### 42-T-2500

Dynamic Condylar Standard Implant Set (Steel Plates in Mat. DIN ISO 5832-1)



### 42-T-2510

Tray (empty) for Dynamic Condylar Standard Implant Set (Steel)

#### Listing:

**42-T-2500** Dynamic Condylar Standard Implant Set (Steel Plates in Mat. DIN ISO 5832-1)

**42-T-2510** Tray for Dynamic Condylar Standard Implant Set (Steel)

Cat. No.:	Dynamic Lag Screw, Thread Length 22 mm	Pcs.:
<b>42-02-0400</b>	Length 50 mm, Steel	2
<b>42-02-0402</b>	Length 55 mm, Steel	2
<b>42-02-0404</b>	Length 60 mm, Steel	2
<b>42-02-0406</b>	Length 65 mm, Steel	2
<b>42-02-0408</b>	Length 70 mm, Steel	2
<b>42-02-0410</b>	Length 75 mm, Steel	1
<b>42-02-0440</b>	Compression Screw	3
Cat. No.:	Dynamic Condylar Screw Plates 95°	Pcs.:
<b>42-02-0450</b>	6 Holes, Length 100 mm, Steel	2
<b>42-02-0452</b>	8 Holes, Length 130 mm, Steel	2
<b>42-02-0454</b>	10 Holes, Length 163 mm, Steel	2
<b>42-02-0456</b>	12 Holes, Length 198 mm, Steel	1
<b>42-02-0458</b>	14 Holes, Length 225 mm, Steel	1
<b>42-02-0460</b>	16 Holes, Length 260 mm, Steel	1




**70-7110-10**  
Sterilization Container (recommended)

## Broad and Narrow Plates Set in Graphic Case

**42-T-0330**

Broad and Narrow Plates Set (Steel Plates in Mat. DIN ISO 5832-1)


**42-T-0331**

Tray (empty) for Broad and Narrow Plates Set

**Listing:**
**42-T-0330** Broad and Narrow Plates Set (Steel Plates in Mat. DIN ISO 5832-1)

**42-T-0331** Tray for Broad and Narrow Plates Set (Steel)

Cat.No.:	Semi-tubular Plates	Pcs.
42-03-0674	4.5 / 4 Hole, Length 71 mm Steel	2
42-03-0676	4.5, 5 / 5 Hole, Length 87 mm Steel	2
42-03-0678	4.5 / 6 Hole, Length 103 mm Steel	2
42-03-0680	4.5 / 7 Hole, Length 119 mm Steel	1

Narrow Dynamic Compression Plates

42-03-0604	4.5 / 4 Hole, Length 71 mm Steel	2
42-03-0606	4.5 / 5 Hole, Length 87 mm Steel	2
42-03-0608	4.5 / 6 Hole, Length 103 mm Steel	4
42-03-0610	4.5 / 7 Hole, Length 119 mm Steel	2
42-03-0612	4.5 / 8 Hole, Length 135 mm Steel	2
42-03-0614	4.5 / 9 Hole, Length 151 mm Steel	2
42-03-0616	4.5 / 10 Hole, Length 167 mm Steel	1
42-03-0620	4.5 / 12 Hole, Length 199 mm Steel	1

Broad Dynamic Compression Plates

42-03-0640	4.5 / 6 Hole, Length 103 mm Steel	2
42-03-0642	4.5 / 7 Hole, Length 119 mm Steel	2
42-03-0644	4.5 / 8 Hole, Length 135 mm Steel	2
42-03-0646	4.5 / 9 Hole, Length 151 mm Steel	1
42-03-0648	4.5 / 10 Hole, Length 167 mm Steel	1
42-03-0652	4.5 / 12 Hole, Length 199 mm Steel	1
42-03-0656	4.5 / 14 Hole, Length 231 mm Steel	1

Spoon Plates

42-03-0692	4.5 / 5 Hole, Length 100 mm Steel	1
42-03-0694	4.5 / 6 Hole, Length 120 mm Steel	1

T-Plates

42-03-0702	4.5 / 4 Hole, Length 84 mm Steel	1
42-03-0706	4.5 / 6 Hole, Length 116 mm Steel	1
42-03-0708	4.5 / 8 Hole, Length 148 mm Steel	1

T-Buttress Plate

42-03-0714	4.5 / 4 Hole, Length 81 mm Steel	1
------------	----------------------------------	---

L-Buttress Plate

42-03-0720	4.5 / 4 Hole, Length 85 mm left Steel	1
42-03-0722	4.5 / 4 Hole, Length 85 mm right Steel	1

**70-7110-10**  
Sterilization Container (recommended)


## Steel Screw Set in Graphic Case



**42-T-0336**  
Screw Set (Steel Screws in Mat. DIN ISO 5832-1)

**42-T-0337**  
Tray (empty) for Screw Set

**Listing:**

**42-T-0336** Screw Plates Set (Steel Screws in Mat. DIN ISO 5832-1)

**42-T-0337** Tray for Screw Set

<u>Cat. No.</u>	<u>Dia 4.5 mm Cortex Screws.</u>	<u>Pcs.</u>
42-03-0200	Length 14 mm	6
42-03-0202	Length 16 mm	6
42-03-0204	Length 18 mm	6
42-03-0206	Length 20 mm	6
42-03-0208	Length 22 mm	6
42-03-0210	Length 24 mm	6
42-03-0212	Length 26 mm	7
42-03-0214	Length 28 mm	7
42-03-0216	Length 30 mm	14
42-03-0218	Length 32 mm	14
42-03-0220	Length 34 mm	14
42-03-0222	Length 36 mm	14
42-03-0224	Length 38 mm	7
42-03-0226	Length 40 mm	14
42-03-0228	Length 42 mm	6
42-03-0230	Length 44 mm	6
42-03-0232	Length 46 mm	6
42-03-0234	Length 48 mm	3
42-03-0238	Length 52 mm	3
42-03-0242	Length 56 mm	3
42-03-0246	Length 60 mm	3
42-03-0248	Length 65 mm	3
42-03-0250	Length 70 mm	3

<u>Cat. No.</u>	<u>Dia 6.5 mm/16 Cortex Screws.</u>	<u>Pcs.</u>
42-03-0396	Length 65 mm	2
42-03-0398	Length 70 mm	2
42-03-0400	Length 75 mm	2
42-03-0402	Length 80 mm	2
42-03-0404	Length 85 mm	2
42-03-0406	Length 90 mm	2
42-03-0408	Length 95 mm	2
42-03-0410	Length 100 mm	2
42-03-0412	Length 105 mm	2
42-03-0414	Length 110 mm	2

<u>Cat. No.</u>	<u>Dia 6.5 mm/16 Cortex Screws.</u>	<u>Pcs.</u>
42-03-0436	Length 45 mm	3
42-03-0438	Length 50 mm	3
42-03-0440	Length 55 mm	3
42-03-0442	Length 60 mm	3
42-03-0444	Length 65 mm	3
42-03-0446	Length 70 mm	3
42-03-0448	Length 75 mm	3
42-03-0450	Length 80 mm	3
42-03-0452	Length 85 mm	3
42-03-0454	Length 90 mm	3
42-03-0456	Length 95 mm	3
42-03-0458	Length 100 mm	3
42-03-0460	Length 105 mm	3
42-03-0462	Length 110 mm	3

<u>Cat. No.</u>	<u>Dia 6.5 mm/16 mm Cancellous Screws. fully threaded</u>	<u>Pcs.</u>
42-03-0480	Length 25 mm	4
42-03-0482	Length 30 mm	4
42-03-0484	Length 35 mm	4
42-03-0486	Length 40 mm	4
42-03-0488	Length 45 mm	2
42-03-0490	Length 50 mm	2
42-03-0492	Length 55 mm	2
42-03-0494	Length 60 mm	2
42-03-0496	Length 65 mm	2
42-03-0498	Length 70 mm	2
42-03-0500	Length 75 mm	2
42-03-0502	Length 80 mm	2
42-03-0504	Length 85 mm	2
42-03-0506	Length 90 mm	2
42-03-0508	Length 95 mm	2
42-03-0510	Length 100 mm	2
42-03-0512	Length 105 mm	2
42-03-0514	Length 110 mm	2

<u>Cat. No.</u>	<u>Dia 4.5 mm Malleolar Screws.</u>	<u>Pcs.</u>
42-03-0352	Length 25/12 mm	3
42-03-0354	Length 30/15 mm	3
42-03-0356	Length 35/18 mm	3
42-03-0358	Length 40/20 mm	3
42-03-0360	Length 45/22 mm	3
42-03-0362	Length 50/24 mm	3
42-03-0364	Length 55/26 mm	3
42-03-0366	Length 60/28 mm	3
42-03-0368	Length 65/30 mm	3
42-03-0370	Length 70/32 mm	3

<u>Cat. No.</u>	<u>Dia 6.5 mm/16 mm Cancellous Screws.</u>	<u>Pcs.</u>
42-03-0382	Length 30 mm	2
42-03-0384	Length 35 mm	2
42-03-0386	Length 40 mm	2
42-03-0388	Length 45 mm	2
42-03-0390	Length 50 mm	2
42-03-0392	Length 55 mm	2
42-03-0394	Length 60 mm	2

<u>Cat. No.</u>	<u>Washer for Dia. 6.5 mm Cancellous Screws</u>	<u>Pcs.</u>
42-03-0522		6

42-04-0184	Forceps	1
------------	---------	---



**70-7310-13**  
Sterilization Container (recommended)

## Titanium Screw Set in Graphic Case



**42-T-0336T**  
Screw Set (Titanium Screws in  
Mat. DIN ISO 5832-3)

**42-T-0337**  
Tray (empty) for  
Screw Set

### Listing:

**42-T-0336T** Screw Plates Set (Titanium Screws in Mat. DIN ISO 5832-3)

**42-T-0337** Tray for Screw Set

Cat. No.	Dia 4.5 mm Cortex Screws.	Pcs.
42-03-0200TI	Length 14 mm	6
42-03-0202TI	Length 16 mm	6
42-03-0204TI	Length 18 mm	6
42-03-0206TI	Length 20 mm	6
42-03-0208TI	Length 22 mm	6
42-03-0210TI	Length 24 mm	6
42-03-0212TI	Length 26 mm	7
42-03-0214TI	Length 28 mm	7
42-03-0216TI	Length 30 mm	14
42-03-0218TI	Length 32 mm	14
42-03-0220TI	Length 34 mm	14
42-03-0222TI	Length 36 mm	14
42-03-0224TI	Length 38 mm	7
42-03-0226TI	Length 40 mm	14
42-03-0228TI	Length 42 mm	6
42-03-0230TI	Length 44 mm	6
42-03-0232TI	Length 46 mm	6
42-03-0234TI	Length 48 mm	3
42-03-0238TI	Length 52 mm	3
42-03-0242TI	Length 56 mm	3
42-03-0246TI	Length 60 mm	3
42-03-0248TI	Length 65 mm	3
42-03-0250TI	Length 70 mm	3

Cat. No.	Dia 4.5 mm Malleolar Screws.	Pcs.
42-03-0352TI	Length 25/12 mm	3
42-03-0354TI	Length 30/15 mm	3
42-03-0356TI	Length 35/18 mm	3
42-03-0358TI	Length 40/20 mm	3
42-03-0360TI	Length 45/22 mm	3
42-03-0362TI	Length 50/24 mm	3
42-03-0364TI	Length 55/26 mm	3
42-03-0366TI	Length 60/28 mm	3
42-03-0368TI	Length 65/30 mm	3
42-03-0370TI	Length 70/32 mm	3

Cat. No.	Dia 6.5 mm/16 mm Cancellous Screws.	Pcs.
42-03-0382TI	Length 30 mm	2
42-03-0384TI	Length 35 mm	2
42-03-0386TI	Length 40 mm	2
42-03-0388TI	Length 45 mm	2
42-03-0390TI	Length 50 mm	2
42-03-0392TI	Length 55 mm	2
42-03-0394TI	Length 60 mm	2

Cat. No.	Dia 6.5 mm/16 Cortex Screws.	Pcs.
42-03-0396TI	Length 65 mm	2
42-03-0398TI	Length 70 mm	2
42-03-0400TI	Length 75 mm	2
42-03-0402TI	Length 80 mm	2
42-03-0404TI	Length 85 mm	2
42-03-0406TI	Length 90 mm	2
42-03-0408TI	Length 95 mm	2
42-03-0410TI	Length 100 mm	2
42-03-0412TI	Length 105 mm	2
42-03-0414TI	Length 110 mm	2

Cat. No.	Dia 6.5 mm/16 Cortex Screws.	Pcs.
42-03-0436TI	Length 45 mm	3
42-03-0438TI	Length 50 mm	3
42-03-0440TI	Length 55 mm	3
42-03-0442TI	Length 60 mm	3
42-03-0444TI	Length 65 mm	3
42-03-0446TI	Length 70 mm	3
42-03-0448TI	Length 75 mm	3
42-03-0450TI	Length 80 mm	3
42-03-0452TI	Length 85 mm	3
42-03-0454TI	Length 90 mm	3
42-03-0456TI	Length 95 mm	3
42-03-0458TI	Length 100 mm	3
42-03-0460TI	Length 105 mm	3
42-03-0462TI	Length 110 mm	3

Cat. No.	Dia 6.5 mm/16 mm Cancellous Screws. fully threaded	Pcs.
42-03-0480TI	Length 25 mm	4
42-03-0482TI	Length 30 mm	4
42-03-0484TI	Length 35 mm	4
42-03-0486TI	Length 40 mm	4
42-03-0488TI	Length 45 mm	2
42-03-0490TI	Length 50 mm	2
42-03-0492TI	Length 55 mm	2
42-03-0494TI	Length 60 mm	2
42-03-0496TI	Length 65 mm	2
42-03-0498TI	Length 70 mm	2
42-03-0500TI	Length 75 mm	2
42-03-0502TI	Length 80 mm	2
42-03-0504TI	Length 85 mm	2
42-03-0506TI	Length 90 mm	2
42-03-0508TI	Length 95 mm	2
42-03-0510TI	Length 100 mm	2
42-03-0512TI	Length 105 mm	2
42-03-0514TI	Length 110 mm	2

Cat. No.		Pcs.
42-03-0522TI	Washer for Dia. 6.5 mm Cancellous Screws	6

42-04-0184	Forceps	1
------------	---------	---





**70-7310-13**  
Sterilization Container (recommended)

## Basic Instrument Set for Screws and Plates in Graphic Case

### 42-T-0340

Basic Instrument Set for Screws and Plates (Steel Plates in Mat. DIN ISO 5832-1)



### 42-T-0341

Tray (empty) for Basic Instrument Set (Steel)

#### Listing:

<b>42-T-0340</b>	Basic Instrument Set for Screws and Plates (Steel Plates in Mat. DIN ISO 5832-1)
<b>42-T-0341</b>	Tray for Basic Instrument Set (Steel)

<u>Cat. No.:</u>		<u>Pcs.:</u>
<b>42-14-1232</b>	Drill Bit , $\varnothing$ 3.2 mm , with end to fit quick coupling	3
<b>42-14-1245</b>	Drill Bit , $\varnothing$ 4.5 mm , with end to fit quick coupling	2
<b>42-09-0104</b>	Malleolar Countersink , 3.2 $\varnothing$ mm tip , for Malleolar Screws	1
<b>42-09-0092</b>	Tap Handle for $\varnothing$ 4.5 mm and $\varnothing$ 6.5 mm Taps and Countersink	1
<b>42-09-0078</b>	Tap , $\varnothing$ 4.5 mm	2
<b>42-09-0079</b>	Tap , $\varnothing$ 6.5 mm	1
<b>42-09-0256</b>	Tap Sleeve , $\varnothing$ 3.5 mm ( also for use as $\varnothing$ 3.2 mm Drill Guide )	1
<b>42-09-0280</b>	Insert Drill Sleeve , $\varnothing$ 4.5 mm / $\varnothing$ 3.2 mm , 58 mm Length	1
<b>42-09-0260</b>	Tap Sleeve , $\varnothing$ 4.5 mm ( also for use as $\varnothing$ 4.5 mm Drill Guide )	1
<b>42-09-0144</b>	Hexagonal Screwdriver , width across flats 3.5 mm	1
<b>42-09-0176</b>	Hexagonal Screwdriver for quick coupling , width across flats 3.5 mm	1
<b>42-09-0114</b>	Depth Gauge for $\varnothing$ 4.5 mm and $\varnothing$ 6.5 mm Screws	1
<b>42-09-0116</b>	Sharp Hook	1
<b>42-09-0262</b>	Drill Sleeve for Plates , 40 mm Length , for Round Hole Plates	1
<b>42-09-0440</b>	Articulated Tension Device with Gauge , 20 mm span	1
<b>42-09-0432</b>	Socket Wrench , width across flats 11 mm	1
<b>42-09-0438</b>	Combination Wrench , width across flats 11 mm	1
<b>42-09-0272</b>	Neutral and Load Drill Guide , $\varnothing$ 4.5 mm	1
<b>42-03-0070</b>	Bending Template , 5 Holes	1
<b>42-03-0071</b>	Bending Template , 7 Holes	1
<b>42-03-0072</b>	Bending Template , 9 Holes	1
<b>42-09-0106</b>	Countersink $\varnothing$ 4.5 mm	1
<b>41-3202-2015</b>	Kirschner Wires $\varnothing$ 2.0 mm / 150 mm ( Steel )	1




**70-7110-10**  
Sterilization Container (recommended)

## Implant Set for Intramedullary Gliding Nails

### 42-T-5000

Implant Set for Intramedullary Gliding Nails (Steel and Titanium Gliding Nail Set)



### 42-T-5001

Tray (empty)

#### Listing:

**42-T-5000** Implant Set for Intramedullary Gliding Nails (Steel and Titanium Gliding Nail Set)

**42-T-5001** Tray (empty)

#### **Steel Intramedullary Gliding Nail Implants ( Stainless Steel DIN ISO 5832-1 )**

<b>42-01-5250</b>	Dia. 1,5 mm x Length 150 mm	4 Pcs.
<b>42-01-5251</b>	Dia. 1,5 mm x Length 400 mm	4 Pcs.
<b>42-01-5252</b>	Dia. 2,0 mm x Length 200 mm	4 Pcs.
<b>42-01-5253</b>	Dia. 2,0 mm x Length 400 mm	4 Pcs.
<b>42-01-5254</b>	Dia. 2,5 mm x Length 250 mm	4 Pcs.
<b>42-01-5255</b>	Dia. 2,5 mm x Length 400 mm	4 Pcs.
<b>42-01-5256</b>	Dia. 3,0 mm x Length 300 mm	4 Pcs.
<b>42-01-5257</b>	Dia. 3,0 mm x Length 400 mm	4 Pcs.
<b>42-01-5258</b>	Dia. 3,5 mm x Length 350 mm	4 Pcs.
<b>42-01-5259</b>	Dia. 3,5 mm x Length 400 mm	4 Pcs.

#### **Titanium Intramedullary Gliding Nail Implants ( Titanium DIN ISO 5832-3 )**

<b>42-01-5263</b>	Dia. 2,0 mm x Length 450 mm Colour Coding Green	4 Pcs.
<b>42-01-5264</b>	Dia. 2,5 mm x Length 450 mm Colour Coding Pink	4 Pcs.
<b>42-01-5265</b>	Dia. 3,0 mm x Length 450 mm Colour Coding Gold	4 Pcs.
<b>42-01-5266</b>	Dia. 3,5 mm x Length 450 mm Colour Coding Blue	4 Pcs.
<b>42-01-5267</b>	Dia. 4,0 mm x Length 450 mm Colour Coding Violet	4 Pcs.
<b>42-01-5268</b>	Dia. 4,5 mm x Length 450 mm Colour Coding Gray	4 Pcs.
<b>42-01-5269</b>	Dia. 5,0 mm x Length 450 mm Colour Coding Silver	4 Pcs.

### 70-7110-10


Sterilization Container (recommended)



## Small Fragment Instrument and Steel Implant Set in Graphic Case

**42-T-0415**

Small Fragment Instrument and Steel Implant Set (Steel Screws and Plates in Mat. DIN ISO 5832-1)



**42-T-0416**

Tray (empty) for Small Fragment Instrument and Steel Implant Set



**70-7110-13**

Sterilization Container (recommended)

## Small Fragment Instrument and Steel Implant Set in Graphic Case

**Listing:**
**42-T-0415** Small Fragment Instrument and Steel Implant Set (Steel Screws and Plates in Mat. DIN ISO 5832-1)

**42-T-0416** Tray for Small Fragment Instrument and Steel Implant Set

<u>Cat.No.</u>	<u>Instruments for Dia. 3.5mm / Dia. 4.0mm Screws</u>	<u>Pcs.:</u>	<u>Cat.No. :</u>	<u>Dia. 4.0 mm Cancellous Screws short Thread</u>	<u>Pcs.:</u>
42-14-1225	Drill Bit Dia. 2.5 mm with quick coupling end	2	42-04-0102	Length 12 mm / 5 mm, Steel	6
42-14-1235	Drill Bit Dia. 3.5 mm with quick coupling end	2	42-04-0104	Length 14 mm / 5 mm, Steel	6
42-09-0102	Small Countersink with quick coupling end	1	42-04-0106	Length 16 mm / 6 mm, Steel	6
42-09-0074	Cortical Tap Dia. 3.5 mm with quick coupling end	2	42-04-0108	Length 18 mm / 7 mm, Steel	6
42-09-0076	Cancellous Tap Dia. 4.0 mm with quick coupling end	2	42-04-0110	Length 20 mm / 8 mm, Steel	6
42-09-0094	Handle with quick coupling	1	42-04-0112	Length 22 mm / 9 mm, Steel	6
42-09-0252	Drill Guide and Drill Sleeve for Drill Bit Dia. 2.0 mm	1	42-04-0114	Length 24 mm / 10 mm, Steel	6
42-09-0254	Tap Sleeve Dia.3.5 mm, Sleeve for Drill Bit Dia.2.5mm	1	42-04-0116	Length 26 mm / 12 mm, Steel	6
42-09-0278	Insert Sleeve Dia.3.5mm / Dia. 2.5mm	1	42-04-0118	Length 28 mm / 14 mm, Steel	6
42-09-0172	Small Hexagonal Screwdriver,insert for Socket 2.5 mm	1	42-04-0120	Length 30 mm / 14 mm, Steel	6
42-09-0174	Small Hexagonal Screwdriver Shaft, Hex. 2.5mm	1	42-04-0122	Length 35 mm / 14 mm, Steel	6
42-09-0112	Small Depth Gauge	1	42-04-0124	Length 40 mm / 14 mm, Steel	6
42-09-0270	Neutral Drill Guide Dia. 3.5mm,used for Drill Bit 2.5mm	1	42-04-0126	Length 45 mm / 15 mm, Steel	6
			42-04-0128	Length 50 mm / 15 mm, Steel	6
<u>Cat.No.:</u>	<u>General Instruments</u>	<u>Pcs.:</u>			
42-09-0116	Sharp Hook	1	42-04-0176	Small Washer for Dia. 3.5 mm / 4.0 mm Screws	6
42-09-0402	Bending Iron for Small Plates	1			
42-09-0404	Bending Iron for Small Plates	1	42-04-0240	Small T-Plate right angle 3/3 Holes Steel	1
42-03-0075	Bending Template for Small Plates 7 Holes	1	42-04-0244	Small T-Plate right angle 4/4 Holes Steel	1
42-03-0076	Bending Template for Small Plates 9 Holes	1	42-04-0250	Small T-Plate oblique angle 3/3 Holes Steel	1
42-09-0580	Self-centring Bone Holding Forceps	1	42-04-0252	Small T-Plate oblique angle 5/3 Holes Steel	1
41-1307	Reduction Forceps with points	1			
42-09-0600	Reduction Forceps for Small Fragments	1	42-04-0220	One Third Tubular Plate 2 Holes, Steel	1
41-0105	Small Hohmann Retractor 6 mm	1	42-04-0222	One Third Tubular Plate 3 Holes, Steel	1
41-0106	Small Hohmann Retractor 8 mm	1	42-04-0224	One Third Tubular Plate 4 Holes, Steel	1
42-09-0654	Hohmann Retractor with broad shank	1	42-04-0226	One Third Tubular Plate 5 Holes, Steel	1
41-0383	Periosteal Elevator, width 6 mm small, straight edge	1	42-04-0228	One Third Tubular Plate 6 Holes, Steel	1
42-04-0184	Screw Forceps self holding	1	42-04-0230	One Third Tubular Plate 7 Holes, Steel	1
			42-04-0232	One Third Tubular Plate 8 Holes, Steel	1
<u>Cat.No.:</u>	<u>Dia. 3.5 mm Cortex Screws</u>	<u>Pcs.:</u>			
42-04-0002	Length 10 mm , Steel	8	42-04-0256	Cloverleaf Plate, Shaft 3 Holes, Steel	1
42-04-0004	Length 12 mm , Steel	8	42-04-0258	Cloverleaf Plate, Shaft 4 Holes, Steel	1
42-04-0006	Length 14 mm , Steel	8			
42-04-0008	Length 16 mm , Steel	8	42-04-0208	Small Fragment Plate 6 Holes, Steel	2
42-04-0010	Length 18 mm , Steel	8	42-04-0210	Small Fragment Plate 7 Holes, Steel	2
42-04-0012	Length 20 mm , Steel	8	42-04-0212	Small Fragment Plate 8 Holes, Steel	2
42-04-0014	Length 22 mm , Steel	8	42-04-0214	Small Fragment Plate 10 Holes, Steel	2
42-04-0016	Length 24 mm , Steel	8			
42-04-0018	Length 26 mm , Steel	8	41-3202-1215	Kirschner Wire Dia.1.2 mm x Length 150 mm, Steel	12
42-04-0020	Length 28 mm , Steel	8	41-3202-1615	Kirschner Wire Dia.1.6 mm x Length 150 mm, Steel	12
42-04-0024	Length 32 mm , Steel	8	41-3202-2015	Kirschner Wire Dia.2.0 mm x Length 150 mm, Steel	12
42-04-0026	Length 36 mm , Steel	8			
42-04-0028	Length 40 mm , Steel	8			
42-04-0030	Length 45 mm , Steel	8			
42-04-0032	Length 50 mm , Steel	8			
<u>Cat.No.:</u>	<u>Dia. 4.0 mm Cancellous Screws full Thread</u>	<u>Pcs.:</u>			
42-04-0144	Length 14 mm, Steel	6			
42-04-0146	Length 16 mm, Steel	6			
42-04-0148	Length 18 mm, Steel	6			
42-04-0150	Length 20 mm, Steel	6			
42-04-0152	Length 22 mm, Steel	6			
42-04-0154	Length 24 mm, Steel	6			
42-04-0156	Length 26 mm, Steel	6			
42-04-0158	Length 28 mm, Steel	6			
42-04-0160	Length 30 mm, Steel	6			
42-04-0162	Length 32 mm, Steel	6			
42-04-0164	Length 35 mm, Steel	6			
42-04-0166	Length 40 mm, Steel	6			
42-04-0168	Length 45 mm, Steel	6			
42-04-0170	Length 50 mm, Steel	6			
42-04-0172	Length 55 mm, Steel	4			
42-04-0174	Length 60 mm, Steel	4			

## Small Fragment Instrument and Titanium Implant Set in Graphic Case

### 42-T-0415T

Small Fragment Instrument and Titanium Implant Set (Titanium Screws and Plates in Mat. DIN ISO 5832-3)



### 42-T-0416

Tray (empty) for Small Fragment Instrument and Titanium Implant Set



### 70-7110-13

Sterilization Container (recommended)

## Small Fragment Instrument and Titanium Implant Set in Graphic Case

**Listing:**
**42-T-0415T** Small Fragment Instrument and Titanium Implant Set (Titanium Screws and Plates in Mat. DIN ISO 5832-3)

**42-T-0416** Tray for Small Fragment Instrument and Titanium Implant Set

<u>Cat.No.</u>	<u>Instruments for Dia. 3.5mm / Dia. 4.0mm Screws</u>	<u>Pcs.:</u>	<u>Cat.No. :</u>	<u>Dia. 4.0 mm Cancellous Screws short Thread</u>	<u>Pcs.:</u>
42-14-1225	Drill Bit Dia. 2.5 mm with quick coupling end	2	42-04-0102TI	Length 12 mm / 5 mm, Titanium	6
42-14-1235	Drill Bit Dia. 3.5 mm with quick coupling end	2	42-04-0104TI	Length 14 mm / 5 mm, Titanium	6
42-09-0102	Small Countersink with quick coupling end	1	42-04-0106TI	Length 16 mm / 6 mm, Titanium	6
42-09-0074	Cortical Tap Dia. 3.5 mm with quick coupling end	2	42-04-0108TI	Length 18 mm / 7 mm, Titanium	6
42-09-0076	Cancellous Tap Dia. 4.0 mm with quick coupling end	2	42-04-0110TI	Length 20 mm / 8 mm, Titanium	6
42-09-0094	Handle with quick coupling	1	42-04-0112TI	Length 22 mm / 9 mm, Titanium	6
42-09-0252	Drill Guide and Drill Sleeve for Drill Bit Dia. 2.0 mm	1	42-04-0114TI	Length 24 mm / 10 mm, Titanium	6
42-09-0254	Tap Sleeve Dia.3.5 mm, Sleeve for Drill Bit Dia.2.5mm	1	42-04-0116TI	Length 26 mm / 12 mm, Titanium	6
42-09-0278	Insert Sleeve Dia.3.5mm / Dia. 2.5mm	1	42-04-0118TI	Length 28 mm / 14 mm, Titanium	6
42-09-0172	Small Hexagonal Screwdriver,insert for Socket 2.5 mm	1	42-04-0120TI	Length 30 mm / 14 mm, Titanium	6
42-09-0174	Small Hexagonal Screwdriver Shaft, Hex. 2.5mm	1	42-04-0122TI	Length 35 mm / 14 mm, Titanium	6
42-09-0112	Small Depth Gauge	1	42-04-0124TI	Length 40 mm / 14 mm, Titanium	6
42-09-0270	Neutral Drill Guide Dia. 3.5mm,used for Drill Bit 2.5mm	1	42-04-0126TI	Length 45 mm / 15 mm, Titanium	6
			42-04-0128TI	Length 50 mm / 15 mm, Titanium	6
<u>Cat.No.:</u>	<u>General Instruments</u>	<u>Pcs.:</u>			
42-09-0116	Sharp Hook	1	42-04-0176TI	Small Washer for Dia. 3.5 mm / 4.0 mm Screws	6
42-09-0402	Bending Iron for Small Plates	1			
42-09-0404	Bending Iron for Small Plates	1	42-04-0240TI	Small T-Plate right angle 3/3 Holes Titanium	1
42-03-0075	Bending Template for Small Plates 7 Holes	1	42-04-0244TI	Small T-Plate right angle 4/4 Holes Titanium	1
42-03-0076	Bending Template for Small Plates 9 Holes	1	42-04-0250TI	Small T-Plate oblique angle 3/3 Holes Titanium	1
42-09-0580	Self-centring Bone Holding Forceps	1	42-04-0252TI	Small T-Plate oblique angle 5/3 Holes Titanium	1
41-1307	Reduction Forceps with points	1			
42-09-0600	Reduction Forceps for Small Fragments	1	42-04-0220TI	One Third Tubular Plate 2 Holes, Titanium	1
41-0105	Small Hohmann Retractor 6 mm	1	42-04-0222TI	One Third Tubular Plate 3 Holes, Titanium	1
41-0106	Small Hohmann Retractor 8 mm	1	42-04-0224TI	One Third Tubular Plate 4 Holes, Titanium	1
42-09-0654	Hohmann Retractor with broad shank	1	42-04-0226TI	One Third Tubular Plate 5 Holes, Titanium	1
41-0383	Periosteal Elevator, width 6 mm small, straight edge	1	42-04-0228TI	One Third Tubular Plate 6 Holes, Titanium	1
42-04-0184	Screw Forceps self holding	1	42-04-0230TI	One Third Tubular Plate 7 Holes, Titanium	1
			42-04-0232TI	One Third Tubular Plate 8 Holes, Titanium	1
<u>Cat.No.:</u>	<u>Dia. 3.5 mm Cortex Screws</u>	<u>Pcs.:</u>			
42-04-0002TI	Length 10 mm , Titanium	8	42-04-0256TI	Cloverleaf Plate, Shaft 3 Holes, Titanium	1
42-04-0004TI	Length 12 mm , Titanium	8	42-04-0258TI	Cloverleaf Plate, Shaft 4 Holes, Titanium	1
42-04-0006TI	Length 14 mm , Titanium	8			
42-04-0008TI	Length 16 mm , Titanium	8	42-04-0208TI	Small Fragment Plate 6 Holes, Titanium	2
42-04-0010TI	Length 18 mm , Titanium	8	42-04-0210TI	Small Fragment Plate 7 Holes, Titanium	2
42-04-0012TI	Length 20 mm , Titanium	8	42-04-0212TI	Small Fragment Plate 8 Holes, Titanium	2
42-04-0014TI	Length 22 mm , Titanium	8	42-04-0214TI	Small Fragment Plate 10 Holes, Titanium	2
42-04-0016TI	Length 24 mm , Titanium	8			
42-04-0018TI	Length 26 mm , Titanium	8			
42-04-0020TI	Length 28 mm , Titanium	8	41-3202-1215	Kirschner Wire Dia.1.2 mm x Length 150 mm, Tit.	12
42-04-0024TI	Length 32 mm , Titanium	8	41-3202-1615	Kirschner Wire Dia.1.6 mm x Length 150 mm, Tit.	12
42-04-0026TI	Length 36 mm , Titanium	8	41-3202-2015	Kirschner Wire Dia.2.0 mm x Length 150 mm, Tit.	12
42-04-0028TI	Length 40 mm , Titanium	8			
42-04-0030TI	Length 45 mm , Titanium	8			
42-04-0032TI	Length 50 mm , Titanium	8			
<u>Cat.No.:</u>	<u>Dia. 4.0 mm Cancellous Screws full Thread</u>	<u>Pcs.:</u>			
42-04-0144TI	Length 14 mm, Titanium	6			
42-04-0146TI	Length 16 mm, Titanium	6			
42-04-0148TI	Length 18 mm, Titanium	6			
42-04-0150TI	Length 20 mm, Titanium	6			
42-04-0152TI	Length 22 mm, Titanium	6			
42-04-0154TI	Length 24 mm, Titanium	6			
42-04-0156TI	Length 26 mm, Titanium	6			
42-04-0158TI	Length 28 mm, Titanium	6			
42-04-0160TI	Length 30 mm, Titanium	6			
42-04-0162TI	Length 32 mm, Titanium	6			
42-04-0164TI	Length 35 mm, Titanium	6			
42-04-0166TI	Length 40 mm, Titanium	6			
42-04-0168TI	Length 45 mm, Titanium	6			
42-04-0170TI	Length 50 mm, Titanium	6			
42-04-0172TI	Length 55 mm, Titanium	4			
42-04-0174TI	Length 60 mm, Titanium	4			

## Mini Instrument and Steel Implant Set in Graphic Case

**42-T-0420**

Mini Instrument and Steel Implant Set (Steel Screws and Plates in Mat. DIN ISO 5832-1)



**42-T-0421**

Tray (empty) for Mini Instrument and Implant Set (Steel)



**70-7110-10**

Sterilization Container (recommended)

## Mini Instrument and Steel Implant Set in Graphic Case

**Listing:**
**42-T-0420** Mini Instrument and Steel Implant Set (Steel Screws and Plates in Mat. DIN ISO 5832-1)

**42-T-0421** Tray for Mini Instrument and Implant Set (Steel)

Cat.No	Instruments for Dia. 2.0mm / Dia. 2.7mm Screws	Pcs.:	Cat.No.:	Dia. 1.5 mm Cortex Screws	Pcs.:
42-14-1220	Drill Bit Dia. 2.0 mm with quick coupling	2	42-04-0400	Length 6 mm , Steel	4
42-14-1227	Drill Bit Dia. 2.7 mm with quick coupling	2	42-04-0402	Length 7 mm , Steel	4
42-09-0102	Small Countersink with end to fit quick coupling Handle	1	42-04-0404	Length 8 mm , Steel	4
42-09-0073	Tap Dia. 2.7 mm with end to fit quick coupling Handle	2	42-04-0406	Length 9 mm , Steel	4
42-09-0094	Handle with quick coupling	1	42-04-0408	Length 10 mm , Steel	4
42-09-0252	Drill Guide and Drill Sleeve for Dia. 2.0 mm Drill Bit	1	42-04-0410	Length 11 mm , Steel	4
42-09-0256	Tap Sleeve Dia. 3.5 mm	1	42-04-0412	Length 12 mm , Steel	4
42-09-0172	Small Hex- Screwdriver 2.5mm with holding Sleeve	1	42-04-0414	Length 14 mm , Steel	4
42-09-0112	Depth Gauge ( only for 2.7 mm Screws )	1	42-04-0416	Length 16 mm , Steel	4

Cat.No.:	Instruments f. Dia.1.5mm/Dia.2.0mm Mini Screws	Pcs.:	Cat.No.:	Dia. 2.0 mm Cortex Screws	Pcs.:
42-14-1211	Drill Bit Dia. 1.1 mm with quick coupling	2	42-04-0430	Length 6 mm , Steel	6
42-14-1215	Drill Bit Dia. 1.5 mm with quick coupling	2	42-04-0432	Length 8 mm , Steel	6
42-09-0100	Mini Countersink to fit Handle with mini quick coupling	1	42-04-0434	Length 10 mm , Steel	6
42-09-0071	Tap Dia. 1.5 mm to fit Handle with mini quick coupling	2	42-04-0436	Length 12 mm , Steel	6
42-09-0072	Tap Dia. 2.0 mm to fit Handle with mini quick coupling	2	42-04-0438	Length 14 mm , Steel	6
42-09-0096	Handle with mini quick coupling	1	42-04-0440	Length 16 mm , Steel	6
42-09-0250	Mini Drill Sleeve for Dia. 1.1 and Dia. 1.5 mm Drill Bits	1	42-04-0442	Length 18 mm , Steel	6
42-09-0110	Mini Depth Gauge ( only for 2.0 mm Screws )	1	42-04-0444	Length 20 mm , Steel	6

Cat.No.:	General Instruments	Pcs.:	Cat.No.:	Plates for Dia. 1.5 mm / 2.0 mm Cortex Screws	Pcs.:
42-09-0116	Sharp Hook	1	42-04-0590	Mini Plate 2.0 straight , 3 Hole / Length 17 mm , Steel	2
42-09-0170	Hexagonal Screw Driver 1.5 mm with holding Sleeve	1	42-04-0592	Mini Plate 2.0 straight , 4 Hole / Length 23 mm , Steel	2
42-09-0400	Bending Iron for Finger Plate	1	42-04-0594	Mini Plate 2.0 straight , 5 Hole / Length 29 mm , Steel	2
41-1309	Holding Forceps for small Plates	1	42-04-0600	Mini L - Plate 2.0 oblique left 2/2 Hole , Steel	2
41-1307	Reduction Forceps with points	1	42-04-0602	Mini L - Plate 2.0 oblique right 2/2 Hole , Steel	2
41-0105	Small Retractor, 6 mm	1	42-04-0604	Mini T - Plate 2.0 , 2/2 Hole , Steel	4
41-0106	Small Retractor, 8 mm	1			
42-09-0654	Retractor for small fragments with broad shank , short	2			
41-0385	Periosteal Elevator, small with straight edge 3 mm	1	42-04-0184	Screwholding forceps	1

Cat.No.:	Dia. 2.7 mm Cortex Screws	Pcs.:
42-04-0460	Length 6 mm , Steel	6
42-04-0462	Length 8 mm , Steel	6
42-04-0464	Length 10 mm , Steel	6
42-04-0466	Length 12 mm , Steel	6
42-04-0468	Length 14 mm , Steel	6
42-04-0470	Length 16 mm , Steel	6
42-04-0472	Length 18 mm , Steel	6
42-04-0474	Length 20 mm , Steel	6
42-04-0476	Length 22 mm , Steel	6
42-04-0478	Length 24 mm , Steel	6

**42-04-0176** Small Washer for Dia. 2.7 mm Screws , Steel 6

Cat.No.:	Plates for Dia. 2.7 mm Cortex Screws	Pcs.:
42-04-0540	Quarter Tubular Plate Length 25 mm / 3 Hole , Steel	2
42-04-0542	Quarter Tubular Plate Length 33 mm / 4 Hole , Steel	2
42-04-0544	Quarter Tubular Plate Length 41 mm / 5 Hole , Steel	2
42-04-0546	Quarter Tubular Plate Length 49 mm / 6 Hole , Steel	2
42-04-0548	Quarter Tubular Plate Length 57 mm / 7 Hole , Steel	2
42-04-0550	Quarter Tubular Plate Length 65 mm / 8 Hole , Steel	2
42-04-0552	L - Plate , oblique left 3/2 Hole , Steel	2
42-04-0554	L - Plate , oblique right 3/2 Hole , Steel	2
42-04-0556	T - Plate , Shaft 3 Hole / Head 2 Hole , Steel	4

Cat.No.:	Kirschner Wires with Trocar	Pcs.:
41-3202-0807	Dia. 0.8 mm / Length 70 mm , Steel	10
41-3202-1007	Dia. 1.0 mm / Length 70 mm , Steel	10
41-3202-1015	Dia. 1.0 mm / Length 150 mm , Steel	10
41-3202-1215	Dia. 1.2 mm / Length 150 mm , Steel	10
41-3202-1615	Dia. 1.6 mm / Length 150 mm , Steel	10

## Mini Instrument and Titanium Implant Set in Graphic Case

### 42-T-0420T

Mini Instrument and Titanium Implant Set (Titanium Screws and Plates in Mat. DIN ISO 5832-3)



### 42-T-0421

Tray (empty) for Mini Instrument and Implant Set (Steel)



### 70-7110-10

Sterilization Container (recommended)

## Mini Instrument and Titanium Implant Set in Graphic Case

**Listing:**
**42-T-0420T** Mini Instrument and Titanium Implant Set (Titanium Screws and Plates in Mat. DIN ISO 5832-3)

**42-T-0421** Tray for Mini Instrument and Implant Set (Steel)

<u>Cat.No</u>	<u>Instruments for Dia. 2.0mm / Dia. 2.7mm Screws</u>	<u>Pcs.:</u>	<u>Cat.No.:</u>	<u>Dia. 1.5 mm Cortex Screws</u>	<u>Pcs.:</u>
<b>42-14-1220</b>	Drill Bit Dia. 2.0 mm with quick coupling	2	<b>42-04-0400TI</b>	Length 6 mm , Titanium	4
<b>42-14-1227</b>	Drill Bit Dia. 2.7 mm with quick coupling	2	<b>42-04-0402TI</b>	Length 7 mm , Titanium	4
<b>42-09-0102</b>	Small Countersink with end to fit quick coupling Handle	1	<b>42-04-0404TI</b>	Length 8 mm , Titanium	4
<b>42-09-0073</b>	Tap Dia. 2.7 mm with end to fit quick coupling Handle	2	<b>42-04-0406TI</b>	Length 9 mm , Titanium	4
<b>42-09-0094</b>	Handle with quick coupling	1	<b>42-04-0408TI</b>	Length 10 mm , Titanium	4
<b>42-09-0252</b>	Drill Guide and Drill Sleeve for Dia. 2.0 mm Drill Bit	1	<b>42-04-0410TI</b>	Length 11 mm , Titanium	4
<b>42-09-0256</b>	Tap Sleeve Dia. 3.5 mm	1	<b>42-04-0412TI</b>	Length 12 mm , Titanium	4
<b>42-09-0172</b>	Small Hex- Screwdriver 2.5mm with holding Sleeve	1	<b>42-04-0414TI</b>	Length 14 mm , Titanium	4
<b>42-09-0112</b>	Depth Gauge ( only for 2.7 mm Screws )	1	<b>42-04-0416TI</b>	Length 16 mm , Titanium	4

<u>Cat.No.:</u>	<u>Instruments for Dia.1.5mm/Dia.2.0mm Mini Screws</u>	<u>Pcs.:</u>	<u>Cat.No.:</u>	<u>Dia. 2.0 mm Cortex Screws</u>	<u>Pcs.:</u>
<b>42-14-1211</b>	Drill Bit Dia. 1.1 mm with quick coupling	2	<b>42-04-0430TI</b>	Length 6 mm , Titanium	6
<b>42-14-1215</b>	Drill Bit Dia. 1.5 mm with quick coupling	2	<b>42-04-0432TI</b>	Length 8 mm , Titanium	6
<b>42-09-0100</b>	Mini Countersink to fit Handle with mini quick coupling	1	<b>42-04-0434TI</b>	Length 10 mm , Titanium	6
<b>42-09-0071</b>	Tap Dia. 1.5 mm to fit Handle with mini quick coupling	2	<b>42-04-0436TI</b>	Length 12 mm , Titanium	6
<b>42-09-0072</b>	Tap Dia. 2.0 mm to fit Handle with mini quick coupling	2	<b>42-04-0438TI</b>	Length 14 mm , Titanium	6
<b>42-09-0096</b>	Handle with mini quick coupling	1	<b>42-04-0440TI</b>	Length 16 mm , Titanium	6
<b>42-09-0250</b>	Mini Drill Sleeve for Dia. 1.1 and Dia. 1.5 mm Drill Bits	1	<b>42-04-0442TI</b>	Length 18 mm , Titanium	6
<b>42-09-0110</b>	Mini Depth Gauge ( only for 2.0 mm Screws )	1	<b>42-04-0444TI</b>	Length 20 mm , Titanium	6

<u>Cat.No.:</u>	<u>General Instruments</u>	<u>Pcs.:</u>	<u>Cat.No.:</u>	<u>Plates for Dia. 1.5 mm / 2.0 mm Cortex Screws</u>	<u>Pcs.:</u>
<b>42-09-0116</b>	Sharp Hook	1	<b>42-04-0590TI</b>	Mini Plate 2.0 straight , 3 Hole / Length 17 mm , Tit.	2
<b>42-09-0170</b>	Hexagonal Screw Driver 1.5 mm with holding Sleeve	1	<b>42-04-0592TI</b>	Mini Plate 2.0 straight , 4 Hole / Length 23 mm , Tit.	2
<b>42-09-0400</b>	Bending Iron for Finger Plate	1	<b>42-04-0594TI</b>	Mini Plate 2.0 straight , 5 Hole / Length 29 mm , Tit.	2
<b>41-1309</b>	Holding Forceps for small Plates	1	<b>42-04-0600TI</b>	Mini L - Plate 2.0 oblique left 2/2 Hole , Titanium	2
<b>41-1307</b>	Reduction Forceps with points	1	<b>42-04-0602TI</b>	Mini L - Plate 2.0 oblique right 2/2 Hole , Titanium	2
<b>41-0105</b>	Small Retractor, 6 mm	1	<b>42-04-0604TI</b>	Mini T – Plate 2.0 , 2/2 Hole , Titanium	4
<b>41-0106</b>	Small Retractor, 8 mm	1			
<b>42-09-0654</b>	Retractor for small fragments with broad shank , short	2			
<b>41-0385</b>	Periosteal Elevator, small with straight edge 3 mm	1	<b>42-04-0184</b>	Screwholding forceps	1

<u>Cat.No.:</u>	<u>Dia. 2.7 mm Cortex Screws</u>	<u>Pcs.:</u>
<b>42-04-0460TI</b>	Length 6 mm , Titanium	6
<b>42-04-0462TI</b>	Length 8 mm , Titanium	6
<b>42-04-0464TI</b>	Length 10 mm , Titanium	6
<b>42-04-0466TI</b>	Length 12 mm , Titanium	6
<b>42-04-0468TI</b>	Length 14 mm , Titanium	6
<b>42-04-0470TI</b>	Length 16 mm , Titanium	6
<b>42-04-0472TI</b>	Length 18 mm , Titanium	6
<b>42-04-0474TI</b>	Length 20 mm , Titanium	6
<b>42-04-0476TI</b>	Length 22 mm , Titanium	6
<b>42-04-0478TI</b>	Length 24 mm , Titanium	6

**42-04-0176TI** Small Washer for Dia. 2.7 mm Screws , Titanium 6



<u>Cat.No.:</u>	<u>Plates for Dia. 2.7 mm Cortex Screws</u>	<u>Pcs.:</u>
<b>42-04-0540TI</b>	Quarter Tubular Plate Length 25 mm / 3 Hole, Titanium	2
<b>42-04-0542TI</b>	Quarter Tubular Plate Length 33 mm / 4 Hole, Titanium	2
<b>42-04-0544TI</b>	Quarter Tubular Plate Length 41 mm / 5 Hole, Titanium	2
<b>42-04-0546TI</b>	Quarter Tubular Plate Length 49 mm / 6 Hole, Titanium	2
<b>42-04-0548TI</b>	Quarter Tubular Plate Length 57 mm / 7 Hole, Titanium	2
<b>42-04-0550TI</b>	Quarter Tubular Plate Length 65 mm / 8 Hole, Titanium	2
<b>42-04-0552TI</b>	L - Plate , oblique left 3/2 Hole , Titanium	2
<b>42-04-0554TI</b>	L - Plate , oblique right 3/2 Hole , Titanium	2
<b>42-04-0556TI</b>	T - Plate , Shaft 3 Hole / Head 2 Hole , Titanium	4

<u>Cat.No.:</u>	<u>Kirschner Wires with Trocar</u>	<u>Pcs.:</u>
<b>41-3202-0807</b>	Dia. 0.8 mm / Length 70 mm , Steel	10
<b>41-3202-1007</b>	Dia. 1.0 mm / Length 70 mm , Steel	10
<b>41-3202-1015</b>	Dia. 1.0 mm / Length 150 mm , Steel	10
<b>41-3202-1215</b>	Dia. 1.2 mm / Length 150 mm , Steel	10
<b>41-3202-1615</b>	Dia. 1.6 mm / Length 150 mm , Steel	10

## Standard Set External Fixateur (Tubular System) in Graphic Case

**42-T-8000**

Graphic Case Standard Set External Fixateur (Tubular System)



**42-T-8001**  
Tray (empty)

**Listing:**

- 42-T-8000**      Graphic Case Standard Set External Fixateur (Tubular System)
- 42-T-8001**      Tray (empty)

<u>Cat.No.</u>		<u>Pcs.:</u>
42-09-0023	Drill Bit Ø 3.5 mm , extra long	2
42-09-0438	Combination Wrench with across flats 11 mm	1
42-08-0068	Simple T-Handle for Steinmann Pins and Schanz Screws Ø5.0 mm	1
42-08-0070	Socket Wrench with across flats 11mm 180mm long for threadet conical Bolts and External Fixateur	1
42-08-0072	Hexagonal Allen Key withacross flats 3.5 mm	1
42-08-0092	Trocar with point Ø 3.5 mm short	1
42-08-0061	Drill Sleeve Ø 5.0 mm / Ø 3.5 mm short	1
42-08-0065	Drill Sleeve Ø 6.0 mm / Ø 5.0 mm short	1
42-08-0062	Trocar with point for Aiming Device Ø 3.5 mm long	1
42-08-0060	Drill Sleeve Ø 5.0 mm / Ø 3.5 mm long	1
42-08-0064	Drill Sleeve Ø 6.0 mm / Ø 5.0 mm long	1

<u>Cat.No.</u>	<u>Fixateur Components</u>	<u>Pcs.:</u>
42-08-0080	Open Compressor	2
42-08-0100	Tube Rod Ø 11 mm , 100 mm Length ( Steel )	1
42-08-0102	Tube Rod Ø 11 mm , 125 mm Length ( Steel )	1
42-08-0104	Tube Rod Ø 11 mm , 150 mm Lenght ( Steel )	2
42-08-0106	Tube Rod Ø 11 mm , 200 mm Length ( Steel )	2
42-08-0108	Tube Rod Ø 11 mm , 250 mm Lenght ( Steel )	2
42-08-0110	Tube Rod Ø 11 mm , 300 mm Length ( Steel )	2
42-08-0112	Tube Rod Ø 11 mm , 350 mm Length ( Steel )	2
42-08-0114	Tube Rod Ø 11 mm , 400 mm Length ( Steel )	2
42-08-0170	Plug for Tube Rods	10
42-08-0188	Open adjustable Clamp	12
42-08-0190	Ø 11 mm Tube-to Tube Clamp	2

**Schanz Screws Ø 5.0 mm , Thread - Length 50 mm ( with Trocar Point and Triangular Shaft )  
( Steel Implants DIN ISO 5832-1 )**

<u>Cat.No.</u>		<u>Pcs.:</u>
42-08-1264	Ø 5.0 mm , Length 125 mm	6
42-08-1266	Ø 5.0 mm , Length 150 mm	6
42-08-1268	Ø 5.0 mm , Length 170 mm	6
42-08-1270	Ø 5.0 mm , Length 200 mm	6



**70-7110-10**  
Sterilization Container (recommended)

## Content

Standard Instrument Set for Dynamic and Condylar Screws and Plates in Graphic Case	42.16.2
Dynamic Hip Screw Standard Implant Set in Graphic Case	42.16.3
Dynamic Condylar Standard Implant Set in Graphic Case	42.16.4
Broad and Narrow Plates Set in Graphic Case	42.16.5
Steel Screw Set in Graphic Case	42.16.6
Titanium Screw Set in Graphic Case	42.16.7
Basic Instrument Set for Screws and Plates in Graphic Case	42.16.8
Implant Set for Intramedullary Gliding Nails	42.16.9
Small Fragment Instrument and Steel Implant Set in Graphic Case	42.16.10
Small Fragment Instrument and Titanium Implant Set in Graphic Case	42.16.12
Mini Instrument and Steel Implant Set in Graphic Case	42.16.14
Mini Instrument and Titanium Implant Set in Graphic Case	42.16.16
Standard Set External Fixateur (Tubular System) in Graphic Case	42.16.18



SCHREIBER GmbH  
D-78567 Fridingen, Germany  
Unterer Damm 15  
P.O. Box 67  
Tel.: +49 7463 307  
Fax.: +49 7463 8870  
E-Mail: [info@schreibergmbh.de](mailto:info@schreibergmbh.de)  
Internet: [www.schreibergmbh.de](http://www.schreibergmbh.de)