





**Sterilcontainer
Sterilizing Containers
Contenedores de Esterilización
Container di Sterilizzazione
Conteneurs de Stérilisation**





Benefits:

- Long lasting sterilization.
- Maximum flexibility in the preparation of surgical sets and sterileable material.
- Protects surgical instruments from damage.
- Standard size allows easy stacking of containers
- Provides safe transportation of surgical instruments.
- Risk free return to storage against contamination.
- Cost-effective.

Features:

All dimensions comply with DIN 58952 standard, thus guaranteeing an economical rate of utilization of the sterilizers.

Container lids and bottoms are made of anodized aluminum-alloy, therefore, due to the thermal conductivity of the material an optimum drying phase is achieved.

All standard size containers can be stacked, thus space and economic storage problems are eliminated.

Wire baskets can also be stacked to prevent valuable instruments from being damaged during preparation of the sets.

Removal containers; ideal for safe return of dirty instruments and textiles without any risk of contamination during transportation after operation.

Wrapping drapes provide a safe sterile medium in the containers and wire baskets; and guarantee the sterilized contents against contamination when they are taken out for operation.

Ensures an economic and long lasting sterilization.

Sterilization:

We recommend following steam sterilization procedures for our containers with perforated lids and non-perforated bottoms, as specified in DIN 58946/6 standard.

For fractional flow procedure:

A few jets of saturated steam of $p_e \geq 1$ bar is generally sufficient to remove the air out of sterilizing chamber.

For fractional vacuum procedure:
















- Repeated evacuation up to a pressure of < 130 mbar in exchange with steam inflow up to a pressure above or below an atmospheric pressure.
- Maintaining of steam inflow until the working overpressure is reached.

For steam injection procedure (Vacuum procedure):

- Single evacuation up to a pressure of < 70 mbar with simultaneous low steam inflow.
- Maintaining of steam inflow until the working pressure is reached.

For further assistance regarding all other gravitation and prevacuum procedures as specified in DIN 58953/9 standard, please refer to the table on next page.

Suitability of our sterile containers for steam sterilization procedures ¹⁾ as per DIN 58953 part 9 (issue Mai 1987)

Sterilization container	Packing		Steam sterilization procedure acc. to DIN 58946, part I				
	Description	Symbolic representation and packing system	Gravitation procedure	Pre-vacuum procedure	Fractional vacuum procedure	Fractional flow procedure	Steam injection procedure
with filters (StB-F)	Textiles in sterilization container with filter in lid and bottom acc. to DIN 58952 part I		+	+	+	+	+
			○	○	+	+	+
	Textiles in sterilization container with filter in lid, bottom not perforated acc. to DIN 58952 part I		-	-	+	+	+
			-	-	+	+	○
	Textiles in sterilization container with filter in bottom, lid not perforated acc. to DIN 58952 part I		○	+	+	+	+
			-	○	+	+	+
	Instruments in wire basket acc. to DIN 58952 part III in sterilization container with filter in lid and bottom acc. to DIN 58952 part I		+ ³⁾	+	+	+	+
	Instruments in wire basket acc. to DIN 58952 part III, in sterilization container with filter in lid and bottom not perforated acc. to DIN 58952 part I		-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
	Instruments in wire basket acc. to DIN 58952 part III, in sterilization container with filter in bottom acc. to DIN 58952 part I		○	+	+	+	+
	with Valves (StB-V)	Textiles in sterilization container with valves in lid and bottom acc. to DIN 58952 part I		-	○	+	+
			-	-	+	+	○
Textiles in sterilization container with valves in lid no valves on the bottom			-	-	+	+	-
			-	-	+	+	-
Instruments in wire basket in sterilization container with valve in lid and bottom			-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
Instruments in wire basket in sterilization container with valve in lid not perforated			-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾

Explanation of symbols: + applicable - not applicable ○ limited application

¹⁾ In order to avoid mistakes with regards to service, operation and/or when loading the sterilizer, please refer to table 1 for correct or wrong combination of sterilization containers for standardized steam sterilization procedures.

²⁾ Suitable are only sterilization containers guaranteeing evaporation or discharge of condensate.

³⁾ Gravitation procedures having no drying device do not guarantee a sufficient after-drying effect.

basic model



Basic model containers are produced with perforated or non-perforated bottom. Filters are not protected against contamination from liquids, dust and other mechanical effects. Lids are available in 6 different colours.



General front view



Each container contains two identification labels; one for surgical department which the container belongs to; one for the name of operation. They are easily slid into their built-in slots and firmly secured by a spring enforcement.



Indicator labels are filled by an authorized surgeon before sterilisation process and inserted into their spring-enforced slots.




The security seal is used to keep the container lock safe from unauthorized opening, and automatically breaks up only when locking flap is lifted.


safe model




Safe containers are fitted with a safety lid protecting filters. They are supplied either with filters or without filters. Safety lids are supplied in 6 different colours.



Locking flaps (1) are movable 180° on a pivot and can be secured in any position. Silicon sealings (2) are provided around the lid, in the center and periphery of the round filter retainer, completely recessed.



The filter retainer is locked by pressing down with a thumb in the direction of (1) and open when pressed with a thumb in the direction of (2).



To lift off container lid: push lock (4) in the direction of arrow and hold lid latch, pull in the direction of arrow (5). For mounting safety lid, simply place it on container lid and press down.

1/1 size containers




1/1 size containers


basic model

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



580 x 280 x 100 mm



-  70-7110-10
-  70-7120-10
-  70-7130-10
-  70-7140-10
-  70-7150-10
-  70-7160-10

-  70-7111-10
-  70-7121-10
-  70-7131-10
-  70-7141-10
-  70-7151-10
-  70-7161-10

580 x 280 x 135 mm





-  70-7110-13
-  70-7120-13
-  70-7130-13
-  70-7140-13
-  70-7150-13
-  70-7160-13

-  70-7111-13
-  70-7121-13
-  70-7131-13
-  70-7141-13
-  70-7151-13
-  70-7161-13

580 x 280 x 150 mm



-  70-7110-15
-  70-7120-15
-  70-7130-15
-  70-7140-15
-  70-7150-15
-  70-7160-15

-  70-7111-15
-  70-7121-15
-  70-7131-15
-  70-7141-15
-  70-7151-15
-  70-7161-15

580 x 280 x 200 mm



-  70-7110-20
-  70-7120-20
-  70-7130-20
-  70-7140-20
-  70-7150-20
-  70-7160-20

-  70-7111-20
-  70-7121-20
-  70-7131-20
-  70-7141-20
-  70-7151-20
-  70-7161-20

580 x 280 x 260 mm



-  70-7110-26
-  70-7120-26
-  70-7130-26
-  70-7140-26
-  70-7150-26
-  70-7160-26


-  70-7111-26
-  70-7121-26
-  70-7131-26
-  70-7141-26
-  70-7151-26
-  70-7161-26

1/1 size containers


safe model

Types and Dimensions

Bottom Non Perforated





Bottom Perforated



580 x 280 x 100 mm






-  70-7112-10
-  70-7122-10
-  70-7132-10
-  70-7142-10
-  70-7152-10
-  70-7162-10

-  70-7113-10
-  70-7123-10
-  70-7133-10
-  70-7143-10
-  70-7153-10
-  70-7163-10

580 x 280 x 135 mm







-  70-7112-13
-  70-7122-13
-  70-7132-13
-  70-7142-13
-  70-7152-13
-  70-7162-13

-  70-7113-13
-  70-7123-13
-  70-7133-13
-  70-7143-13
-  70-7153-13
-  70-7163-13

580 x 280 x 150 mm



-  70-7112-15
-  70-7122-15
-  70-7132-15
-  70-7142-15
-  70-7152-15
-  70-7162-15

-  70-7113-15
-  70-7123-15
-  70-7133-15
-  70-7143-15
-  70-7153-15
-  70-7163-15

580 x 280 x 200 mm





-  70-7112-20
-  70-7122-20
-  70-7132-20
-  70-7142-20
-  70-7152-20
-  70-7162-20

-  70-7113-20
-  70-7123-20
-  70-7133-20
-  70-7143-20
-  70-7153-20
-  70-7163-20

580 x 280 x 260 mm



-  70-7112-26
-  70-7122-26
-  70-7132-26
-  70-7142-26
-  70-7152-26
-  70-7162-26


-  70-7113-26
-  70-7123-26
-  70-7133-26
-  70-7143-26
-  70-7153-26
-  70-7163-26

1/1 size containers

bottoms and lids


Types and Dimensions

Bottom Non Perforated



580 x 280 x 100 mm	70-7105-10
580 x 280 x 135 mm	70-7105-13
580 x 280 x 150 mm	70-7105-15
580 x 280 x 200 mm	70-7105-20
580 x 280 x 260 mm	70-7105-26

Bottom Perforated



580 x 280 x 100 mm	70-7106-10
580 x 280 x 135 mm	70-7106-13
580 x 280 x 150 mm	70-7106-15
580 x 280 x 200 mm	70-7106-20
580 x 280 x 260 mm	70-7106-26

Types and Dimensions

Lid Perforated



580 x 280 mm




70-7100-01
70-7100-02
70-7100-03
70-7100-04
70-7100-05
70-7100-06

Safety Lid



580 x 280 mm




70-7103-01
70-7103-02
70-7103-03
70-7103-04
70-7103-05
70-7103-06

accessories


Special Accessories

Wire Baskets



485 x 255 x 30 mm **70-7190-03**
 485 x 255 x 50 mm **70-7190-05**
 485 x 255 x 70 mm **70-7190-07**
 485 x 255 x 100 mm **70-7190-10**


Wire Baskets



540 x 255 x 30 mm **70-7191-03**
 540 x 255 x 50 mm **70-7191-05**
 540 x 255 x 70 mm **70-7191-07**
 540 x 255 x 100 mm **70-7191-10**

Special Accessories

Wire Base




545 x 255 x 30 mm **70-7193-00**

Drape Retainer




545 x 255 mm **70-7195-00**

Drape



130 x 60 cm **70-7095-07**
 140 x 100 cm **70-7095-08**
 185 x 45 cm **70-7095-09**
 185 x 60 cm **70-7095-10**

Silicon Mat



520 x 230 mm **70-7852-30**
 440 x 230 mm **70-7852-20**

3/4 size containers




3/4 size containers


basic model

Types and Dimensions

Bottom Non Perforated










Bottom Perforated



465 x 280 x 100 mm



-  70-7210-10
-  70-7220-10
-  70-7230-10
-  70-7240-10
-  70-7250-10
-  70-7260-10

-  70-7211-10
-  70-7221-10
-  70-7231-10
-  70-7241-10
-  70-7251-10
-  70-7261-10

465 x 280 x 135 mm






-  70-7210-13
-  70-7220-13
-  70-7230-13
-  70-7240-13
-  70-7250-13
-  70-7260-13

-  70-7211-13
-  70-7221-13
-  70-7231-13
-  70-7241-13
-  70-7251-13
-  70-7261-13

465 x 280 x 150 mm



-  70-7210-15
-  70-7220-15
-  70-7230-15
-  70-7240-15
-  70-7250-15
-  70-7260-15


-  70-7211-15
-  70-7221-15
-  70-7231-15
-  70-7241-15
-  70-7251-15
-  70-7261-15

3/4 size containers


safe model

Types and Dimensions

Bottom Non Perforated









Bottom Perforated



465 x 280 x 100 mm





-  70-7212-10
-  70-7222-10
-  70-7232-10
-  70-7242-10
-  70-7252-10
-  70-7262-10

-  70-7213-10
-  70-7223-10
-  70-7233-10
-  70-7243-10
-  70-7253-10
-  70-7263-10

465 x 280 x 135 mm








-  70-7212-13
-  70-7222-13
-  70-7232-13
-  70-7242-13
-  70-7252-13
-  70-7262-13

-  70-7213-13
-  70-7223-13
-  70-7233-13
-  70-7243-13
-  70-7253-13
-  70-7263-13

465 x 280 x 150 mm



-  70-7212-15
-  70-7222-15
-  70-7232-15
-  70-7242-15
-  70-7252-15
-  70-7262-15


-  70-7213-15
-  70-7223-15
-  70-7233-15
-  70-7243-15
-  70-7253-15
-  70-7263-15

3/4 size containers


bottoms and lids

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



465 x 280 x 100 mm **70-7205-10**
 465 x 280 x 135 mm **70-7205-13**
 465 x 280 x 150 mm **70-7205-15**



465 x 280 x 100 mm **70-7206-10**
 465 x 280 x 135 mm **70-7206-13**
 465 x 280 x 150 mm **70-7206-15**

Types and Dimensions

Lid Perforated



Safety Lid



70-7200-01
70-7200-02
70-7200-03
70-7200-04
70-7200-05
70-7200-06




70-7203-01
70-7203-02
70-7203-03
70-7203-04
70-7203-05
70-7203-06

3/4 size containers accessories


Special Accessories

Wire Baskets



405 x 255 x 30 mm **70-7290-03**
405 x 255 x 50 mm **70-7290-05**
405 x 255 x 70 mm **70-7290-07**
405 x 255 x 100 mm **70-7290-10**

Wire Base



405 x 255 x 30 mm **70-7293-03**


Special Accessories

Drape Retainer




405 x 255 mm **70-7295-00**

Drape



140 x 100 cm **70-7095-08**

Silicon Mat



440 x 230 mm **70-7852-20**
380 x 230 mm **70-7852-10**

1/2 size containers




1/2 size containers


basic model

Types and Dimensions

Bottom Non Perforated




Bottom Perforated



285 x 280 x 100 mm



-  70-7310-10
-  70-7320-10
-  70-7330-10
-  70-7340-10
-  70-7350-10
-  70-7360-10

-  70-7311-10
-  70-7321-10
-  70-7331-10
-  70-7341-10
-  70-7351-10
-  70-7361-10

285 x 280 x 135 mm




-  70-7310-13
-  70-7320-13
-  70-7330-13
-  70-7340-13
-  70-7350-13
-  70-7360-13

-  70-7311-13
-  70-7321-13
-  70-7331-13
-  70-7341-13
-  70-7351-13
-  70-7361-13

285 x 280 x 150 mm



-  70-7310-15
-  70-7320-15
-  70-7330-15
-  70-7340-15
-  70-7350-15
-  70-7360-15

-  70-7311-15
-  70-7321-15
-  70-7331-15
-  70-7341-15
-  70-7351-15
-  70-7361-15

285 x 280 x 200 mm




-  70-7310-20
-  70-7320-20
-  70-7330-20
-  70-7340-20
-  70-7350-20
-  70-7360-20

-  70-7311-20
-  70-7321-20
-  70-7331-20
-  70-7341-20
-  70-7351-20
-  70-7361-20

285 x 280 x 260 mm



-  70-7310-26
-  70-7320-26
-  70-7330-26
-  70-7340-26
-  70-7350-26
-  70-7360-26


-  70-7311-26
-  70-7321-26
-  70-7331-26
-  70-7341-26
-  70-7351-26
-  70-7361-26

1/2 size containers


safe model

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



285 x 280 x 100 mm



-  70-7312-10
-  70-7322-10
-  70-7332-10
-  70-7342-10
-  70-7352-10
-  70-7362-10

-  70-7313-10
-  70-7323-10
-  70-7333-10
-  70-7343-10
-  70-7353-10
-  70-7363-10

285 x 280 x 135 mm



-  70-7312-13
-  70-7322-13
-  70-7332-13
-  70-7342-13
-  70-7352-13
-  70-7362-13

-  70-7313-13
-  70-7323-13
-  70-7333-13
-  70-7343-13
-  70-7353-13
-  70-7363-13

285 x 280 x 150 mm



-  70-7312-15
-  70-7322-15
-  70-7332-15
-  70-7342-15
-  70-7352-15
-  70-7362-15

-  70-7313-15
-  70-7323-15
-  70-7333-15
-  70-7343-15
-  70-7353-15
-  70-7363-15

285 x 280 x 200 mm




-  70-7312-20
-  70-7322-20
-  70-7332-20
-  70-7342-20
-  70-7352-20
-  70-7362-20

-  70-7313-20
-  70-7323-20
-  70-7333-20
-  70-7343-20
-  70-7353-20
-  70-7363-20

285 x 280 x 260 mm



-  70-7312-26
-  70-7322-26
-  70-7332-26
-  70-7342-26
-  70-7352-26
-  70-7362-26


-  70-7313-26
-  70-7323-26
-  70-7333-26
-  70-7343-26
-  70-7353-26
-  70-7363-26

1/2 size containers

Bottom and Lids


Types and Dimensions

Bottom Non Perforated



285 x 280 x 100 mm	70-7305-10
285 x 280 x 135 mm	70-7305-13
285 x 280 x 150 mm	70-7305-15
285 x 280 x 200 mm	70-7305-20
285 x 280 x 260 mm	70-7305-26

Bottom Perforated



285 x 280 x 100 mm	70-7306-10
285 x 280 x 135 mm	70-7306-13
285 x 280 x 150 mm	70-7306-15
285 x 280 x 200 mm	70-7306-20
285 x 280 x 260 mm	70-7306-26

Types and Dimensions

Lid Perforated



285x 280 mm



70-7300-01
70-7300-02
70-7300-03
70-7300-04
70-7300-05
70-7300-06

Safety Lid



285 x 280 mm




70-7303-01
70-7303-02
70-7303-03
70-7303-04
70-7303-05
70-7303-06

1/2 size containers

accessories


Special Accessories

Wire Baskets



255 x 245 x 30 mm	70-7390-03
255 x 245 x 50 mm	70-7390-05
255 x 245 x 70 mm	70-7390-07
255 x 245 x 100 mm	70-7390-10

Wire Base



258 x 258 x 30 mm **70-7393-00**


Special Accessories

Drape Retainer



255 x 255 mm **70-7395-00**

Drape



90 x 60 cm	70-7095-04
90 x 80 cm	70-7095-05
130 x 45 cm	70-7095-06
130 x 60 cm	70-7095-17

Silicon Mat



245 x 245 mm	70-7852-30
220 x 230 mm	70-7852-00

other containers




Wide body containers


basic model

Types and Dimensions

Bottom Non Perforated







Bottom Perforated



600 x 400 x 120 mm











-  70-7110-62
-  70-7120-62
-  70-7130-62
-  70-7140-62
-  70-7150-62
-  70-7160-62

-  70-7111-62
-  70-7121-62
-  70-7131-62
-  70-7141-62
-  70-7151-62
-  70-7161-62

600 x 400 x 180 mm



-  70-7110-68
-  70-7120-68
-  70-7130-68
-  70-7140-68
-  70-7150-68
-  70-7160-68


-  70-7111-68
-  70-7121-68
-  70-7131-68
-  70-7141-68
-  70-7151-68
-  70-7161-68

Wide body containers


bottoms, lids and accessories

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



Lid Perforated



70-7100-51

70-7100-52

70-7100-53

70-7100-54

70-7100-55

70-7100-56

600 x 400 x 120 mm
600 x 400 x 180 mm


70-7105-62
70-7105-68

600 x 400 x 120 mm
600 x 400 x 180 mm

70-7106-62
70-7106-88


Special Accessories

Wire Baskets



470 x 337 x 51 mm

70-7190-55




533 x 337 x 101 mm

70-7191-60






Types and Dimensions

Bottom Non Perforated








285 x 280 x 55 mm




-  70-7310-05
-  70-7320-05
-  70-7330-05
-  70-7340-05
-  70-7350-05






285x 280 x 85 mm



-  70-7310-08
-  70-7320-08
-  70-7330-08
-  70-7340-08
-  70-7350-08

Bottom Perforated



-  70-7311-08
-  70-7321-08
-  70-7331-08
-  70-7341-08
-  70-7351-08

Bottoms, lids and accessories

Types and Dimensions


Bottom Non Perforated



Bottom Perforated



Lid Perforated



285 x 280 x 55 mm
285 x 280 x 85 mm

70-7305-05
70-7305-08



285 x 280 x 85 mm **70-7306-08**



70-7304-01

70-7304-02


70-7304-03

70-7304-04

70-7304-05

Special Accessories


Wire Baskets



255 x 245 x 30 mm
255 x 245 x 50 mm

70-7390-03
70-7390-05

Wire Base



245 x 255 x 30 mm **70-7393-00**

Silicon Mat



245 x 245 mm
220 x 230 mm


70-7851-00
70-7852-00

Drape Retainer



255 x 255 mm **70-7395-00**

Drape




90 x 60 cm
90 x 80 cm


70-7095-04
70-7095-05

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



300 x 140 x 40 mm






-  70-7410-04
-  70-7420-04
-  70-7430-04
-  70-7440-04
-  70-7450-04

-  70-7411-04
-  70-7421-04
-  70-7431-04
-  70-7441-04
-  70-7451-04

300 x 140 x 70 mm







-  70-7410-07
-  70-7420-07
-  70-7430-07
-  70-7440-07
-  70-7450-07

-  70-7411-07
-  70-7421-07
-  70-7431-07
-  70-7441-07
-  70-7451-07

300 x 140 x 100 mm




-  70-7410-10
-  70-7420-10
-  70-7430-10
-  70-7440-10
-  70-7450-10

-  70-7411-10
-  70-7421-10
-  70-7431-10
-  70-7441-10
-  70-7451-10


bottoms, lids and accessories

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



Lid Perforated



300 x 140 x 40 mm **70-7405-04**
 300 x 140 x 70 mm **70-7405-07**
 300 x 140 x 100 mm **70-7405-10**



300 x 140 x 40 mm **70-7406-04**
 300 x 140 x 70 mm **70-7406-07**
 300 x 140 x 100 mm **70-7406-10**

300 x 140 mm



70-7400-01
70-7400-02
70-7400-03
70-7400-04
70-7400-05

Special Accessories

Silicon Mat




270 x 125 mm **70-7850-00**

Drape Retainer



280 x 125 mm **70-7495-00**

Drape




60 x 40 cm **70-7095-01**

Dental containers


bottoms, lids and accessories

Types and Dimensions

Lid and Bottom Perforated



Bottom Perforated



Lid Perforated



310 x 190 x 40 mm



70-7510-04

310 x 190 x 65 mm




70-7510-06

310 x 190 x 130 mm



70-7510-13



310 x 190 x 40 mm **70-7506-04** 310 x 190 mm **70-7500-01**
 310 x 190 x 65 mm **70-7506-06**
 310 x 190 x 130 mm **70-7506-13**


Special Accessories

Drape Retainer



310 x 190 mm **70-7595-00**

Drape



60 x 40 cm **70-7095-01**
 65 x 65 cm **70-7095-02**
 85 x 85 cm **70-7095-03**

Types and Dimensions

Lid Perforated - Bottom Perforated (Non Filter)



300 x 138 x 65 mm



- 70-7610-06
- 70-7620-06
- 70-7630-06
- 70-7640-06
- 70-7650-06

500 x 155 x 75 mm



- 70-7612-07
- 70-7622-07
- 70-7632-07
- 70-7642-07
- 70-7652-07

Scopy containers

Types and Dimensions

Lid Perforated - Bottom Perforated (Non Filter)



250 x 60 x 50 mm



- 70-7600-01

450 x 70 x 70 mm




- 70-7600-02


Bottom, lids and accessories

Types and Dimensions

Lid And Bottom Non Perforated



Bottom Non Perforated



285 x 280 x 100 mm



285 x 280 x 150 mm



580 x 280 x 150 mm



- 70-7315-10
- 70-7325-10
- 70-7335-10
- 70-7345-10
- 70-7355-10
- 70-7365-10

- 70-7315-15
- 70-7325-15
- 70-7335-15
- 70-7345-15
- 70-7355-15
- 70-7365-15

- 70-7115-15
- 70-7125-15
- 70-7135-15
- 70-7145-15
- 70-7155-15
- 70-7165-15



285 x 280 x 150 mm **70-7305-10**

285 x 280 x 150 mm **70-7305-15**



580 x 280 x 150 mm **70-7105-15**


Lids - Non Perforated



Lids - Non Perforated



580 x 280 mm




- 70-7102-01**
- 70-7102-02**
- 70-7102-03**
- 70-7102-04**
- 70-7102-05**
- 70-7102-06**

285 x 280 mm



- 70-7302-01**
- 70-7302-02**
- 70-7302-03**
- 70-7302-04**
- 70-7302-05**
- 70-7302-06**



70-7880-00

Brush container

285 x 280 x 85 mm

Containers are supplied together with mounting fixtures onto the wall.

For 20 brushes. Brushes are taken out by moving of the handle in both directions right or left which is under the container.



70-7880-00

70-7880-01

70-7880-02

70-7880-03

70-7880-04

Filters, Indicators And Labels

Disposable paper filters



70-7800-00
100 Pcs/Package
For 1/1, 3/4 and Half Size Containers



70-7802-00
100 Pcs/Package
For Mini Containers

70-7804-00
100 Pcs/Package
For Dental Containers



70-7801-00
Filter Dispenser

Textile Filters



70-7800-01
Single
For 1/1, 3/4 and Half Size Containers



70-7802-01
Single
For Mini Containers

70-7804-01
Single
For Dental Containers



70-7809-00
For Paper Filter

70-7809-01
For Textile Filter




70-7810-00
100 Pcs / Package
Paper label with indicator
For 1/1, 3/4 and Half Size Containers



70-7812-00
100 Pcs / Package
Paper label with indicator
For Mini and Dental Containers




70-7819-00
100 Pcs / Package
Plastic Security Seal



- 70-7820-01
- 70-7820-02
- 70-7820-03
- 70-7820-04
- 70-7820-05
- 70-7820-06

Coloured Identification Labels
Engraved to your specifications, max. 16 letters.



70-7829-00
Dispenser for identification labels

Components





70-7830-00
Holding Pin
For Wire Baskets

70-7830-01
Holding Pin
For Partition Sheets



Partition Sheets
Stainless Steel


70-7831-05 50 x 20 mm
70-7831-13 130 x 20 mm
70-7831-22 225 x 20 mm
70-7831-46 460 x 20 mm



70-7840-00
Divider, to be used with 70-7841-00



70-7841-00 232 mm
70-7841-01 471 mm
Sliding Rails with Securing Screws




70-7842-00 232 mm
70-7842-01 471 mm
Limiting Bar with Securing Screws



70-7845-00
Spiral Holding



70-7847-00
Double Spiral Holding



Fixation Clamps



70-7849-00 50 x 20 x 12
70-7849-00 65 x 20 x 20



70-7849-20 Ø 7-12 mm
70-7849-21 Ø 12-16 mm
70-7849-22 Ø 16-26 mm
70-7849-23 Ø 26-36 mm
70-7849-24 Ø 36-45 mm



Carrying



70-7890-00
Trolley
720 x 620 x 800 mm



70-7891-00
Trolley
1030 x 620 x 1170 mm

SCHREIBER GmbH
Unterer Damm 15
D78567 Fridingen
P.O. Box 67
Tel. +49-7463-307
Fax. +49-7463-8870
e-mail info@schreiber GmbH.de
Internet <http://www.schreiber GmbH.de>



Distributed by: